

خەملاندن بۆ ئامپىر بۆ راگە ياندن

ژمارە	ناوى ئامپىر	ژمارەى ئامپىر
1	Hood class 2 include heap filter with UV light size(160cm,LX*90cm,H)	2
2	vortex	2
3	Heater block	1
4	Microtube centrifuge 2 4position	2
5	Micropipet set(Thermoscientefic)	2set
6	T.S.B meter	2
7	Interferon gamma Releas Assay(IGRA-TB-TEST)	
8	Full automated extraction	1
9	Centerfuge(24tube)	20
10	Microscope	20



بايۆلۇجىت/ جنان احمد علي.
دوا به روار بۆ پيشكهش كردنى ئوفەر (۲۰۲۰/۱۱/۶)